

Table 6: Cumulative^a Rates of Inflammation (Safety Population of All Treated Eyes)

Category	Inflammation	Fermented N=211		Control N=210	
		n	%	n	%
Epithelial Edema	None (0)	131	62.1	134	63.8
	Trace (+1)	58	27.5	54	25.7
	Mild (+2)	26	12.3	22	10.5
	Moderate (+3)	7	3.3	7	3.3
	Severe (+4)	1	0.5	1	0.5
Stromal Edema	None (0)	161	76.3	160	76.2
	Trace (+1)	46	21.8	45	21.4
	Mild (+2)	8	3.8	11	5.2
	Moderate (+3)	1	0.5	2	1.0
	Severe (+4)	0	0.0	1	0.5
Cells	None (0)	12	5.7	8	3.8
	Rare (0.5+)	133	63.0	132	62.9
	Trace (+1)	155	73.5	147	70.0
	Mild (+2)	29	13.7	36	17.1
	Moderate (+3)	4	1.9	6	2.9
Flare	None (0)	100	47.4	96	45.7
	Faint/Trace/Mild (1+)	111	52.6	111	52.9
	Moderate (2+)	2	0.9	8	3.8
Anterior Synechiae	None (0)	210	99.5	209	99.5
	Trace (+1)	1	0.5	0	0.0
	Mild (+2)	0	0.0	1	0.5
Posterior Synechiae	None (0)	210	99.5	210	100.0
	Trace (+1)	1	0.5	0	0.0
Fibrin Presence	None (0)	208	98.6	209	99.5
	Trace (+1)	2	0.9	1	0.5
	Mild (+2)	1	0.5	0	0.0

^a Cumulative includes unscheduled visits

Endothelial Cell Count

The primary effectiveness endpoint of the mean percent endothelial cell (ECC) change preoperatively vs. 3 months was achieved in the study as the lower 95% CI of the mean percent difference between OVD groups was within the -5% non-inferiority margin for ECC loss in both the ITT and safety populations (Table 7). As such, the mean percent ECC change in the fermented group was statistically non-inferior to that of the control group.

Table 7: Mean Percent Change Difference in Endothelial Cell Counts

Population	Mean Percent Change ^a	Lower 95% Confidence Limit	Upper 95% Confidence Limit
TT: All Randomized Subjects	1.05 (0.80 SE)	-0.53	2.64
Safety: All Paired-Eye Subjects	1.11 (1.89 SD)	-0.52	2.74

^a Percent Change=(Postop ECC Minus Preop ECC)/Preop ECC with Difference= Fermented OVD Percent Change Minus Control OVD

As shown in Table 8 between preop and 3 months postoperative, the mean percent ECC change was -5.55% (SD 9.99) in the fermented group and -6.66% (SD 10.23) in the control group, for a mean percent change difference between groups of 1.11% (SD 11.89).

Table 8: Change in ECC from Baseline to 3 Months (Safety Population- Paired Eye Subjects)

Variable	Group	N ^a	Mean	Standard Deviation
Preoperative ^b	Fermented	206	2552.96	361.65
	Control	206	2543.75	355.64
	Difference	206	9.21	131.54
3 Months	Fermented	206	2410.82	420.16
	Control	206	2377.14	423.53
	Difference	206	33.68	313.36
Percent Change ^c (Preoperative vs 3 months)	Fermented	206	-5.55	9.99
	Control	206	-6.66	10.23
	Difference	206	1.11 ^d	11.89

^a Two paired-eye subjects missing 3-month ECC data; therefore, N=206.
^b The difference between OVD groups at preoperative visit was not statistically significant (p-value=0.3162)
^c Percent Change=(Postop ECC Minus Preop ECC)/Preop ECC
^d 95% CI: -0.52, 2.74

How supplied

The Healon GV[®] PRO OVD is a sterile, non-pyrogenic, viscoelastic preparation supplied in disposable glass syringes, delivering 0.85 mL sodium hyaluronate (18 mg/mL) dissolved in physiological sodium chloride phosphate buffer (pH 6.8 – 7.6). Each mL of Healon GV[®] PRO OVD contains 18 mg of sodium hyaluronate.

The Healon GV[®] PRO OVD syringes are terminally sterilized and aseptically packaged. A sterile, single-use 27G cannula is enclosed in the 0.85mL blister.

Preparation and storage
 Refrigerated Healon GV[®] PRO OVD should be held at room temperature for approximately 30 minutes before use. Protect from freezing and exposure to light.

For intraocular use.

Store between 2 to 8°C (36 to 46°F).

Definition of symbols

	Caution, see instructions for use
	See instructions for use
	Do not reuse
	Protect from light
	Do not use if the packaging has been opened or damaged
	Protect from freezing
	Temperature limitation 2°C (36°F) to 8°C (46°F)
	Sterilized using steam (solution)
	Sterilized using irradiation (cannula)
	Sterilized using aseptic processing techniques (blister packaging)
	Manufacturer
	Batch code
	Use by (YYYY-MM-DD: year-month-day)
	Catalogue number

References

- Balazs, E.A.: Ultrapure hyaluronic acid and the use thereof. U.S. patent 4,141,973 (1979)
- Fry L.L. & Yee R.W. (1993): Healon GV in extracapsular cataract extraction with intraocular lens implantation. *Cataract Refract. Surg.* 19:409-412.
- Gaskel A. & Haining W. (1991): A double blind randomized multicentre clinical trial of "Healon GV," compared with "Healon" in ECCE with IOL implantation. *Eur J. Implant Ref. Surg.* 3:241.
- Pape, L.G. & Balazs, E.A.: The use of sodium hyaluronate (Healon[®]) in human anterior segment surgery. *Ophthalmol* 87 (1980) p 699-705

Rx only

AMO Uppsala AB
 Rapskatan 7
 Box 6406
 SE-751 36 Uppsala, Sweden

Product of Sweden

Johnson & Johnson Surgical Vision, Inc.
 1700 E. St. Andrew Place
 Santa Ana, CA 92705 USA
 1-877-266-4543

Revision date: 11/2020

HEALON GV and HEALON5 are trademarks of Johnson & Johnson Surgical Vision, Inc. All rights reserved.



Healon GV[®] PRO

Sodium Hyaluronate


NOT MADE WITH NATURAL RUBBER LATEX

Instructions

Sterile opening technique

Tear off the paper covering.

Remove syringe and place on sterile field.




Assembly


Press the vial completely into the holder so that the needle perforates the membrane.

Important


Perforate the membrane before screwing on the plastic rod.




Remove the plastic rod.



Screw the plastic rod into the plunger.
 Connect the cannula.



Check for proper function.



Store at 2 to 8°C (36 to 46°F).

For single use only.



